

Acquisition Method Info

Method Name Joan_silvia_MRM_HILIC_02102019_labeled_5ul.m
 Method Path D:\MassHunter\Methods\Joan_silvia_MRM_HILIC_02102019_labeled_5ul.m
 Method Description Veure volulm i dilució de la mostra

Device List

HiP Sampler
 Binary Pump
 Column Comp.
 QQQ

MS QQQ Mass Spectrometer

Ion Source AJS ESI Tune File D:\MassHunter\Tune\QQQ\G6490A\atunes.TUNE.XML
 Stop Mode No Limit/As Pump Stop Time (min) 1
 Time Filter On Time Filter Width (min) 0.07

Time Segments

Index	Start Time (min)	Scan Type	Ion Mode	Div Valve	Delta EMV	Store
1	0	MRM	ESI+Agilent Jet Stream	To Waste	300	Yes
2	0.35	MRM	ESI+Agilent Jet Stream	To MS	300	Yes

Time Segment 1

Scan Segments

Cpd Name	ISTD?	Prec Ion	MS1 Res	Prod Ion	MS2 Res	Dwell	Frag (V)	CE (V)	Cell Acc (V)	Polarity
Compound 1	No	350	Wide / Unit (6490)	200	Wide / Unit (6490)	10	380	0	5	Positive

Scan Parameters

Data Stg Threshold
 Centroid 0

Source Parameters

Parameter	Value (+)	Value (-)
Gas Temp (°C)	270	270
Gas Flow (l/min)	15	15
Nebulizer (psi)	35	35
SheathGasHeater	400	400
SheathGasFlow	11	11
Capillary (V)	3000	3000
VCharging	1000	1500

Ion Funnel Parameters

Pos High Pressure RF	130	Neg High Pressure RF	110
Pos Low Pressure RF	100	Neg Low Pressure RF	60

Time Segment 2

Acquisition Method Report



Scan Segments

Cpd Name	ISTD?	Prec Ion	MS1 Res	Prod Ion	MS2 Res	Dwell	Frag (V)	CE (V)	Cell Acc (V)	Polarity
Ac_CoA_2	No	819	Wide / Unit (6490)	435	Wide / Unit (6490)	10	380	25	5	Positive
Ac_CoA_3	No	819	Wide / Unit (6490)	305	Wide / Unit (6490)	10	380	36	5	Positive
Ac_CoA_1	No	819	Wide / Unit (6490)	138	Wide / Unit (6490)	10	380	60	5	Positive
Ac_CoA_1_4	No	817	Wide / Unit (6490)	435	Wide / Unit (6490)	10	380	25	5	Positive
Ac_CoA_2_0	No	817	Wide / Unit (6490)	433	Wide / Unit (6490)	10	380	25	5	Positive
Ac_CoA_2_1	No	817	Wide / Unit (6490)	305	Wide / Unit (6490)	10	380	36	5	Positive
Ac_CoA_1_5	No	817	Wide / Unit (6490)	303	Wide / Unit (6490)	10	380	36	5	Positive
Ac_CoA_1_3	No	817	Wide / Unit (6490)	138	Wide / Unit (6490)	10	380	60	5	Positive
Ac_CoA_1_9	No	817	Wide / Unit (6490)	136	Wide / Unit (6490)	10	380	60	5	Positive
Ac_CoA_8	No	815	Wide / Unit (6490)	433	Wide / Unit (6490)	10	380	25	5	Positive
Ac_CoA_9	No	815	Wide / Unit (6490)	303	Wide / Unit (6490)	10	380	36	5	Positive
Ac_CoA_7	No	815	Wide / Unit (6490)	136	Wide / Unit (6490)	10	380	60	5	Positive
Ac_CoA_1_7	No	814	Wide / Unit (6490)	430	Wide / Unit (6490)	10	380	25	5	Positive
Ac_CoA_1_8	No	814	Wide / Unit (6490)	305	Wide / Unit (6490)	10	380	36	5	Positive
Ac_CoA_1_6	No	814	Wide / Unit (6490)	138	Wide / Unit (6490)	10	380	60	5	Positive
Ac_CoA_5	No	812	Wide / Unit (6490)	430	Wide / Unit (6490)	10	380	25	5	Positive
Ac_CoA_1_1	No	812	Wide / Unit (6490)	428	Wide / Unit (6490)	10	380	25	5	Positive
Ac_CoA_1_2	No	812	Wide / Unit (6490)	305	Wide / Unit (6490)	10	380	36	5	Positive
Ac_CoA_6	No	812	Wide / Unit (6490)	303	Wide / Unit (6490)	10	380	36	5	Positive
Ac_CoA_4	No	812	Wide / Unit (6490)	138	Wide / Unit (6490)	10	380	60	5	Positive
Ac_CoA_1_0	No	812	Wide / Unit (6490)	136	Wide / Unit (6490)	10	380	60	5	Positive
Ac_CoA_3_u	No	810	Wide / Unit (6490)	428	Wide / Unit (6490)	10	380	25	5	Positive
Ac_CoA_1_u	No	810	Wide / Unit (6490)	303	Wide / Unit (6490)	10	380	36	5	Positive
Ac_CoA_2_u	No	810	Wide / Unit (6490)	136	Wide / Unit (6490)	10	380	60	5	Positive
SAM_2	No	407	Wide / Unit (6490)	306	Wide / Unit (6490)	10	380	4	5	Positive
SAM_1	No	407	Wide / Unit (6490)	102	Wide / Unit (6490)	10	380	32	5	Positive
SAM_10	No	406	Wide / Unit (6490)	305	Wide / Unit (6490)	10	380	4	5	Positive
SAM_9	No	406	Wide / Unit (6490)	102	Wide / Unit (6490)	10	380	32	5	Positive
SAM_14	No	405	Wide / Unit (6490)	304	Wide / Unit (6490)	10	380	4	5	Positive
SAM_13	No	405	Wide / Unit (6490)	102	Wide / Unit (6490)	10	380	32	5	Positive
SAM_6	No	404	Wide / Unit (6490)	303	Wide / Unit (6490)	10	380	4	5	Positive
SAM_5	No	404	Wide / Unit (6490)	102	Wide / Unit (6490)	10	380	32	5	Positive
SAM_12	No	402	Wide / Unit (6490)	301	Wide / Unit (6490)	10	380	4	5	Positive
SAM_11	No	402	Wide / Unit (6490)	97	Wide / Unit (6490)	10	380	32	5	Positive
SAM_4	No	401	Wide / Unit (6490)	300	Wide / Unit (6490)	10	380	4	5	Positive
SAM_3	No	401	Wide / Unit (6490)	97	Wide / Unit (6490)	10	380	32	5	Positive
SAM_8	No	400	Wide / Unit (6490)	299	Wide / Unit (6490)	10	380	4	5	Positive
SAM_7	No	400	Wide / Unit (6490)	97	Wide / Unit (6490)	10	380	32	5	Positive
SAM_3u	No	399	Wide / Unit (6490)	298	Wide / Unit (6490)	10	380	4	5	Positive
SAM_1u	No	399	Wide / Unit (6490)	250	Wide / Unit (6490)	10	380	12	5	Positive
SAM_2u	No	399	Wide / Unit (6490)	97	Wide / Unit (6490)	10	380	32	5	Positive

Scan Parameters

Data Stg	Threshold
Centroid	0

Source Parameters

Parameter	Value (+)	Value (-)
Gas Temp (°C)	270	270
Gas Flow (l/min)	15	15
Nebulizer (psi)	35	35
SheathGasHeater	400	400
SheathGasFlow	11	11
Capillary (V)	3000	3000
VCharging	1000	1500

Ion Funnel Parameters

Pos High Pressure RF	130	Neg High Pressure RF	110
Pos Low Pressure RF	100	Neg Low Pressure RF	60

Chromatograms

Chrom Type	Label	Offset	Y-Range
TIC	TIC	0	10000000

Instrument Curves

Actual

Name: HiP Sampler

Model: G4226A

Auxiliary

Draw Speed	50.0 µL/min
Eject Speed	100.0 µL/min
Draw Position Offset	0.0 mm
Wait Time After Drawing	2.0 s
Sample Flush Out Factor	5.0
Vial/Well bottom sensing	Yes

Injection

Injection Mode	Standard injection
Injection Volume	5.00 µL

High throughput

Automatic Delay Volume Reduction	No
Overlapped Injection	
Enable Overlapped Injection	No

Valve Switching

Valve Movements	0
Valve Switch Time 1	
Switch Time 1 Enabled	No
Valve Switch Time 2	
Switch Time 2 Enabled	No
Valve Switch Time 3	
Switch Time 3 Enabled	No
Valve Switch Time 4	
Switch Time 4 Enabled	No

Stop Time

Stoptime Mode	As pump/No limit
---------------	------------------

Post Time

Posttime Mode	Off
---------------	-----

Name: Binary Pump Model: G4220A

Flow 0.400 mL/min
Use Solvent Types Yes
Stroke Mode Synchronized
Low Pressure Limit 0.00 bar
High Pressure Limit 900.00 bar
Max. Flow Ramp Up 100.000 mL/min²
Max. Flow Ramp Down 100.000 mL/min²
Expected Mixer No check

Stroke A

Automatic Stroke Calculation A Yes

Stop Time

Stoptime Mode Time set
Stoptime 13.00 min

Post Time

Posttime Mode Off

Solvent Composition

	Channel	Ch. 1 Solv.	Name 1	Ch2 Solv.	Name 2	Selected	Used	Percent
1	A	100.0 % Water V.03		100.0 % Water V.03		Ch. 1	Yes	2.00 %
2	B	100.0 % Acetonitrile V.03		100.0 % Acetonitrile V.03		Ch. 1	Yes	98.00 %

Timetable

	Time	A	B	Flow	Pressure
1	2.00 min	2.00 %	98.00 %	0.400 mL/min	900.00 bar
2	9.00 min	60.00 %	40.00 %	0.400 mL/min	900.00 bar
3	9.50 min	2.00 %	98.00 %	0.400 mL/min	900.00 bar
4	13.00 min	2.00 %	98.00 %	0.400 mL/min	900.00 bar

Name: Column Comp. Model: G1316C

Valve Position Position 1 (Port 1 -> 2)
Ready when front door open Yes

Left Temperature Control

Temperature Control Mode Temperature Set
Temperature 25.0 °C

Enable Analysis Left Temperature

Enable Analysis Left Temperature On Yes
Enable Analysis Left Temperature Value 0.8 °C

Right Temperature Control

Right temperature Control Mode Temperature Set
Right temperature 25.0 °C

Enable Analysis Right Temperature

Enable Analysis Right Temperature On Yes
Enable Analysis Right Temperature Value 0.8 °C

Stop Time

Stoptime Mode As pump/injector

Post Time

Posttime Mode Off